

Enterovirus Ev71 3D Rabbit Monoclonal Antibody (Wb, Elisa) - Ev71 3D

Item Number: CM0004980

Introduction

Enterovirus EV71 3D Rabbit Monoclonal Antibody (mAb) for Western blot (WB) and ELISA applications. Recombinantly expressed immunogen; detects EV71 3D protein with human cross-reactivity. Ideal for virology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	EV71 3D Rabbit mAb
Species	Enterovirus
Immunogen	Recombinant protein of Enterovirus EV71 3D
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Protein Weight	51kDa
Shipping	Ice bag