


Erab/Hsd17B10 Rabbit Monoclonal Antibody For Wb, Elisa - Q99714

Item Number: CM0002201



Introduction

High-quality rabbit monoclonal antibody against human ERAB/HSD17B10 (HSD17B10, gene ID 3028), suitable for WB and ELISA. Cross-reacts with human, mouse, rat samples. Detects a 27kDa protein band.

[Learn More](#)

Parameter	Value
Product Name	ERAB/HSD17B10 Rabbit mAb
Remarks/Alias	ABAD; CAMR; ERAB; HCD2; MHBD; HADH2; MRPP2; MRX17; MRX31; SCHAD; MRXS10; SDR5C1; HSD10MD; 17b-HSD10; DUPXp11.22; ERAB/HSD17B10
Species	Human
Gene ID (Human)	3028
Gene ID	3028
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 162-261 of human ERAB/HSD17B10 (Q99714).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99714
Protein Weight	27kDa
Shipping	Ice bag