


Eral1 Rabbit Polyclonal Antibody For Wb And Elisa - 075616

Item Number: CM0016446



Introduction

High-quality Anti-ERAL1 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat ERAL1, a mitochondrial GTPase critical for 12S rRNA binding and small subunit assembly. Ideal for mitochondrial research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	ERAL1 Rabbit pAb
Remarks/Alias	ERA; CEGA; ERA-W; H-ERA; ERAL1A; ERAL1B; HERA-A; HERA-B; PRLTS6; ERAL1
Species	Human
GeneID (Human)	26284
GeneID	26284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 198-437 of human ERAL1 (NP_005693.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75616
Protein Weight	48kDa
Shipping	Ice bag