


Esra Rabbit Mab For Wb - P03372

Item Number: CM0004061



Introduction

Recombinant ESRα rabbit monoclonal antibody validated for WB, IP, ELISA, ChIP, and Flow Cytometry. Targets human estrogen receptor alpha (ERα) for hormone signaling, cancer, and reproductive research.

[Learn More](#)

Parameter	Value
Product Name	ESRα Rabbit mAb
Remarks/Alias	ER; ESRα
Species	Human
Gene ID (Human)	2099
Immunogen	Synthetic peptide corresponding to amino acids 495-595 of human ESRα (NP_000116.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA, ChIP, FC (intra)
Cross Reactivity	Human
SWISS	P03372
Protein Weight	66kDa
Shipping	Ice bag