


Esrp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Nxxg1

Item Number: CM0014038



Introduction

ESRP1 Rabbit Polyclonal Antibody validated for WB and ELISA in human and mouse. Targets epithelial splicing regulatory protein 1, a splicing factor involved in EMT and epithelial-specific isoform regulation. Suitable for research in deafness and cancer.

[Learn More](#)

Parameter	Value
Product Name	ESRP1 Rabbit pAb
Remarks/Alias	RBM35A; RMB35A; DFNB109; ESRP1
Species	Human
Gene ID	54845
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ESRP1 (NP_060167.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6NXG1
Protein Weight	76kDa
Shipping	Ice bag