


Esrrb Rabbit Polyclonal Antibody For Wb And Elisa - 095718

Item Number: CM0027990



Introduction

ESRRB rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human and mouse ERR2 nuclear receptor. Suitable for stem cell pluripotency, early development, and transcriptional regulation research. 095718.

[Learn More](#)

Parameter	Value
Product Name	ESRRB Rabbit pAb
Remarks/Alias	ERR2; ERRb; ESRL2; NR3B2; DFNB35; ERRbeta2; ERR beta-2; ESRRB
Species	Human
GeneID (Human)	2103
GeneID	2103
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 434-508 of human ESRRB (NP_004443.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95718
Protein Weight	48kDa
Shipping	Ice bag