

Etf1 Rabbit Pab For Wb, If/Icc, Elisa - P62495

Item Number: CM0030520



Introduction

Polyclonal rabbit antibody against human ETF1, validated for WB, IF/ICC, ELISA. Recognizes human, mouse, and rat ETF1 (49 kDa), a translation termination factor. Ideal for studying protein synthesis and nonsense-mediated decay.

[Learn More](#)

Parameter	Value
Product Name	ETF1 Rabbit pAb
Remarks/Alias	ERF; RF1; ERF1; TB3-1; D5S1995; SUP45L1; ETF1
Species	Human
GeneID (Human)	2107
GeneID	2107
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 278-437 of human ETF1 (NP_004721.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62495
Protein Weight	49kDa
Shipping	Ice bag