

Etnk2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nvf9

Item Number: CM0017830

Introduction

High-specificity ETNK2 Rabbit Polyclonal Antibody validated for Western blotting and ELISA. Detects human ethanolamine kinase 2, involved in phospholipid biosynthesis. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	ETNK2 Rabbit pAb
Remarks/Alias	EKI2; HMFT1716; ETNK2
Species	Human
Gene ID (human)	55224
Gene ID	55224
Immunogen	Recombinant protein or recombinant fusion protein containing a sequence corresponding to amino acids 1-386 of human ETNK2 (NP_060678.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NVF9
Protein Weight	45kDa
Shipping	Ice bag