


Ets1 Rabbit Pab - P14921

Item Number: CM0019093



Introduction

Rabbit polyclonal antibody against human ETS1 transcription factor for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with human and mouse. Suitable for gene regulation, immunity, and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ETS1 Rabbit pAb
Remarks/Alias	p54; ETS-1; EWSR2; c-ets-1; ETS1
Species	Human
Gene ID	2113 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ETS1 (NP_005229.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14921
Protein Weight	50kDa
Shipping	Ice bag