

Etv4 Rabbit Polyclonal Antibody For Wb, Elisa - P43268

Item Number: CM0020755



Introduction

ETV4 Rabbit Polyclonal Antibody (P43268) validated for WB and ELISA. Reacts with human, mouse, and rat. Recombinant immunogen covering amino acids 40-246 of human ETV4. Targets a transcriptional activator involved in motor neuron positioning and keratinocyte differentiation.

[Learn More](#)

Parameter	Value
Product Name	ETV4 Rabbit pAb
Remarks/Alias	E1AF; PEA3; E1A-F; PEAS3; ETV4
Species	Human
GeneID (Human)	2118
GeneID	2118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-246 of human ETV4 (NP_001977.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43268
Protein Weight	54kDa
Shipping	Ice bag