

Ev71 3C Protease Rabbit Polyclonal Antibody For Wb And Elisa - Ev71 3C Rabbit Pab

Item Number: CM0021432

Introduction

Polyclonal rabbit antibody targeting the EV71 3C protease (aa 1-183). Validated for Western blot and ELISA with predicted wide cross-reactivity. Suitable for enterovirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	EV71 3C Rabbit pAb
Species	Enterovirus
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human EV71 3C. (B9VUU3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Other (Wide Range Predicted)
Shipping	Ice bag