

Ev71 3C Rabbit Mab - Ev71 3C

Item Number: CM0006569

Introduction

Rabbit monoclonal antibody against enterovirus EV71 3C protease, validated for WB and ELISA. Suitable for enterovirus detection and research, with a predicted molecular weight of 20 kDa.

[Learn More](#)

Parameter	Specification
Product Name	EV71 3C Rabbit mAb
Species	Enterovirus
Immunogen	Recombinant protein of enterovirus EV71 3C
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Enterovirus
Protein Weight	20KDa
Shipping	Ice bag