

Ev71 3D Rabbit Pab - B9Vuu3

Item Number: CM0020366

Introduction

EV71 3D Rabbit pAb is a polyclonal antibody raised in rabbit against Enterovirus EV71 3D protein (aa 1-217, B9VUU3). Suitable for WB and ELISA applications with human cross-reactivity. Ideal for EV71 research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	EV71 3D Rabbit pAb
Species	Enterovirus
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-217 of Enterovirus EV71 3D (B9VUU3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Shipping	Ice bag