

Ev71 3D Rabbit Polyclonal Antibody - Ev71 3D

Item Number: CM0030785

Introduction

Rabbit polyclonal antibody raised against recombinant Enterovirus EV71 3D protein. Applicable in WB and ELISA, with reported human cross-reactivity. Supports EV71 detection and assay development.

[Learn More](#)

Parameter	Value
Product Name	EV71 3D Rabbit pAb
Species	Enterovirus
Immunogen	Recombinant protein of Enterovirus EV71 3D
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Shipping	Ice bag