

Ewsr1 Rabbit Monoclonal Antibody For Wb, Ihc-P - Q01844

Item Number: CM0010300



Introduction

EWSR1 Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Targets human EWSR1 N-terminus; reacts with human, mouse, rat. Ideal for cancer and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	EWSR1 Rabbit mAb
Remarks/Alias	EWS; EWS-FLI1; bK984G1.4; EWSR1
Species	Human
Gene ID (Human)	2130
Gene ID	2130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human EWSR1 (Q01844)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01844
Protein Weight	68kDa
Shipping	Ice bag