

Ewsr1 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q01844

Item Number: CM0016801



Introduction

Rabbit polyclonal antibody against human EWSR1, validated for IF/ICC and ELISA applications. Targets the Ewing sarcoma breakpoint region 1 protein, a transcriptional repressor involved in tumorigenesis.

[Learn More](#)

Parameter	Value
Product Name	EWSR1 Rabbit pAb
Remarks/Alias	EWS; EWS-FLI1; bK984G1.4; EWSR1
Species	Human
GeneID (Human)	2130
GeneID	2130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EWSR1 (NP_001156757.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q01844
Protein Weight	68kDa
Shipping	Ice bag