

Exoc3 Rabbit Polyclonal Antibody For Wb, Elisa - O60645

Item Number: CM0021013

Introduction

High-quality rabbit polyclonal antibody against EXOC3 (Sec6), validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Suitable for exocyst complex and exocytosis research in diagnostic and laboratory settings.

[Learn More](#)

Parameter	Value
Product Name	EXOC3 Rabbit pAb
Remarks/Alias	SEC6; Sec6p; SEC6L1; EXOC3
Species	Human
Gene ID (Human)	11336
Gene ID	11336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 396-745 of human EXOC3 (NP_009208.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60645
Protein Weight	86kDa
Shipping	Ice bag