

Exosc4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npd3

Item Number: CM0029152



Introduction

EXOSC4 Rabbit Polyclonal Antibody validated for WB and ELISA against human and rat samples, targeting the RNA exosome component RRP41 (Q9NPD3). Ideal for RNA processing and degradation studies.

[Learn More](#)

Parameter	Value
Product Name	EXOSC4 Rabbit pAb
Remarks/Alias	SKI6; p12A; RRP41; Ski6p; RRP41A; Rrp41p; hRrp41p; EXOSC4
Species	Human
GeneID (Human)	54512
GeneID	54512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human EXOSC4 (NP_061910.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NPD3
Protein Weight	26kDa
Shipping	Ice bag