


Extl3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 043909

Item Number: CM0030676



Introduction

Rabbit polyclonal antibody against human EXTL3, a glycosyltransferase and REG3 receptor. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for heparan sulfate and REG3 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	EXTL3 Rabbit pAb
Remarks/Alias	RPR; BOTV; REGR; EXTR1; ISDNA; EXTL1L; EXTL3
Species	Human
Gene ID (Human)	2137
Gene ID	2137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-220 of human EXTL3 (NP_001431.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43909
Protein Weight	105kDa
Shipping	Ice bag