


F2R Rabbit Pab - P25116

Item Number: CM0024730



Introduction

High-quality anti-F2R rabbit polyclonal antibody for Western blot and ELISA detection of human, mouse, and rat proteinase-activated receptor 1. Ideal for coagulation research and cardiovascular studies.

[Learn More](#)

Parameter	Value
Product Name	F2R Rabbit pAb
Remarks/Alias	TR; HTR; CF2R; PAR1; PAR-1; F2R
Species	Human
Gene ID (Human)	2149
Gene ID	2149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-110 of human F2R (NP_001983.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25116
Protein Weight	47kDa
Shipping	Ice bag