

F4/80 Rabbit Mab For If, Ihc-P, Elisa - Q61549

Item Number: CM0006997

Introduction

Rabbit monoclonal antibody against mouse F4/80 (Emr1), suitable for immunofluorescence, immunohistochemistry (paraffin), and ELISA. Detects mouse and rat F4/80, a macrophage marker involved in cell adhesion and Treg development. Ideal for immunological research.

[Learn More](#)

Parameter	Value
Product Name	F4/80 Rabbit mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7
Species	Mouse
GeneID (Human)	2015
GeneID	13733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-645 of mouse F4/80 (NP_034260.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag