

Fadd Rabbit Polymab® Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q61160 / Q13158

Item Number: CM0009805



Introduction

High-specificity FADD Rabbit PolymAb® monoclonal antibody for apoptosis research. Validated in WB, IP, IHC-P, ELISA with cross-reactivity to human, mouse, and rat FADD (Gene ID 8772/14082).

[Learn More](#)

Parameter	Value
Product Name	FADD Rabbit PolymAb®
Remarks/Alias	Mort1/FADD
Species	Human/Mouse
GeneID (Human)	8772
GeneID (Mouse)	14082
Immunogen	Recombinant protein (recombinant fusion protein of human/mouse FADD)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q61160/Q13158
Protein Weight	23kDa
Shipping	Ice bag