


Faf1 Rabbit Polyclonal Antibody (Wb/Elisa) - Q9Unn5

Item Number: CM0027020



Introduction

Rabbit polyclonal antibody against human FAF1 (FAS-associated factor 1), suitable for Western blot and ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for apoptosis and DNA replication studies.

[Learn More](#)

Parameter	Value
Product Name	FAF1 Rabbit pAb
Remarks/Alias	hFAF1; CGI-03; HFAF1s; UBXD12; UBXN3A; FAF1
Species	Human
GeneID (Human)	11124
GeneID	11124
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human FAF1 (NP_008982.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNN5
Protein Weight	74kDa
Shipping	Ice bag