

Fam160B2 Rabbit Pab - Q80Yr2

Item Number: CM0019182

Introduction

FAM160B2 Rabbit pAb is a polyclonal antibody for WB, IHC-P, IF/ICC, IP, ELISA. Validated in human, mouse, rat, targeting FHF complex subunit HOOK-interacting protein 2B involved in MAPK/ERK signaling and intestinal barrier function.

[Learn More](#)

Parameter	Value
Product Name	FAM160B2 Rabbit pAb
Remarks/Alias	Rai16; Fam160b2; G430067P06Rik; FAM160B2
Species	Mouse
GeneID (Human)	239170
GeneID	239170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 74-303 of mouse FAM160B2 (NP_919326.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q80YR2
Protein Weight	82kDa
Shipping	Ice bag