


## Fam50A Rabbit Mab (Anti-Fam50A Monoclonal Antibody) For Wb, Ihc-P, Elisa - Q14320

Item Number: CM0007163



### Introduction

Rabbit monoclonal antibody targeting human FAM50A, validated for western blot, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat proteins. Ideal for studying pre-mRNA splicing regulation.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | FAM50A Rabbit mAb                                                                                                                     |
| Remarks/Alias    | XAP5; HXC26; MRXSA; HXC-26; DXS9928E; FAM50A                                                                                          |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 9130                                                                                                                                  |
| GeneID           | 9130                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FAM50A (NP_004690.1) |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | WB, IHC-P, ELISA                                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | Q14320                                                                                                                                |
| Protein Weight   | 40kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |