

Fancd2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Bxw9

Item Number: CM0022952



Introduction

Rabbit polyclonal antibody against FANCD2, validated for Western blot, immunofluorescence, immunoprecipitation, and ELISA. Detects human, mouse, rat FANCD2 (164kDa). Suitable for DNA repair and Fanconi anemia research.

[Learn More](#)

Parameter	Value
Product Name	FANCD2 Rabbit pAb
Remarks/Alias	FA4; FAD; FADC; FAD2; FA-D2; FANCD; FANCD2
Species	Human
Gene ID (Human)	2177
Gene ID	2177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-291 of FANCD2 (NP_001018125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXW9
Protein Weight	164kDa
Shipping	Ice bag