


Faslg Rabbit Pab - P48023

Item Number: CM0014159



Introduction

FasLG rabbit polyclonal antibody for WB, IHC-P, ELISA. Recognizes human and mouse FasLG, a key apoptosis mediator. Ideal for studying immune cell death and immune tolerance.

[Learn More](#)

Parameter	Value
Product Name	FasLG Rabbit pAb
Remarks/Alias	APTL, FASL, CD178, CD95L, ALPS1B, CD95-L, TNFSF6, TNLG1A, APT1LG1, FasLG
Species	Human
Gene ID	356
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 101-200 of human FasLG (NP_000630.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48023
Protein Weight	31kDa
Shipping	Ice bag