

Fast Myosin Skeletal Heavy Chain Rabbit Mab - P12882

Item Number: CM0009205

Introduction

Rabbit monoclonal antibody against human MYH1 (Myosin-1, skeletal muscle heavy chain). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for muscle physiology and hearing research.

[Learn More](#)

Parameter	Value
Product Name	Fast Myosin Skeletal Heavy chain Rabbit mAb
Remarks/Alias	MYHa; HEL71; MYHSA1; MyHC-2x; MyHC-2X/D
Species	Human
GeneID (Human)	4619
GeneID	4619
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human MYH1 (NP_005954.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12882
Protein Weight	223kDa
Shipping	Ice bag