

Fastkd5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7L8L6

Item Number: CM0024532

Introduction

Rabbit polyclonal antibody targeting FASTKD5, suitable for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Ideal for mitochondrial RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	FASTKD5 Rabbit pAb
Remarks/Alias	dj1187M17.5; FASTKD5
Species	Human
GeneID (Human)	60493
GeneID	60493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 495-764 of human FASTKD5 (NP_068598.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L8L6
Protein Weight	87kDa
Shipping	Ice bag