

Fau Rabbit Pab - P35544

Item Number: CM0017575

Introduction

FAU rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications. Detects endogenous FAU/RPS30 protein (~14 kDa) across human, mouse, and rat species. Ideal for apoptosis and ubiquitin-like fusion protein research.

[Learn More](#)

Parameter	Specification
Product Name	FAU Rabbit pAb
Remarks/Alias	S30; FAU1; Fub1; Fubi; asr1; eS30; RPS30; MNSFbeta; FAU
Species	Human
Gene ID	2197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-133 of human FAU (NP_001988.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35544
Protein Weight	14kDa
Shipping	Ice bag