

Fbln1 Rabbit Polyclonal Antibody For Wb, Elisa - P23142

Item Number: CM0023090

Introduction

High-quality rabbit polyclonal antibody against human Fibulin-1 (FBLN1), validated for WB and ELISA. Reacts with human, mouse, and rat FBLN1. Suitable for extracellular matrix and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	FBLN1 Rabbit pAb
Remarks/Alias	FBLN; FIBL1; FBLN1
Species	Human
Gene ID (Human)	2192
Gene ID	2192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 570-670 of human FBLN1 (NP_001987.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23142
Protein Weight	77kDa
Shipping	Ice bag