

Fcgrt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55899

Item Number: CM0025764



Introduction

High-quality rabbit polyclonal antibody against human FCGRT (FcRn), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for IgG recycling and neonatal immunity research.

[Learn More](#)

Parameter	Value
Product Name	FCGRT Rabbit pAb
Remarks/Alias	FCRN; FcγgammaRn; alpha-chain; FCGRT
Species	Human
Gene ID (Human)	2217
Gene ID	2217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-290 of human FcRn (NP_004098.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55899
Protein Weight	40kDa
Shipping	Ice bag