

Fchsd1 Rabbit Pab - Q86Wn1

Item Number: CM0017951

Introduction

Rabbit polyclonal antibody against human FCHSD1. Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for studying actin polymerization and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	FCHSD1 Rabbit pAb
Remarks/Alias	NWK2; FCHSD1
Species	Human
GeneID (Human)	89848
GeneID	89848
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human FCHSD1 (NP_258260.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WN1
Protein Weight	77kDa
Shipping	Ice bag