

Fcrl3 Rabbit Polyclonal Antibody For Wb And Elisa - Q96P31

Item Number: CM0016029

Introduction

Rabbit polyclonal antibody against human FCRL3 (CD307c, MAIA), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for B-cell biology, immunology, and fertilization research.

[Learn More](#)

Parameter	Value
Product Name	FCRL3 Rabbit pAb
Remarks/Alias	MAIA; FCRH3; IFGP3; IRTA3; SPAP2; CD307c; FCRL3
Species	Human
GeneID (human)	115352
GeneID	115352
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-720 of human FCRL3 (NP_443171.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P31
Protein Weight	81kDa
Shipping	Ice bag