


Fes Rabbit Polyclonal Antibody For Wb, Elisa - P07332

Item Number: CM0027649



Introduction

FES Rabbit Polyclonal Antibody targets human FES (93 kDa) with cross-reactivity to mouse. Suitable for WB and ELISA applications in tyrosine kinase signaling, mast cell degranulation, and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	FES Rabbit pAb
Remarks/Alias	FPS; FES
Species	Human
GeneID (Human)	2242
GeneID	2242
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-420 of human FES (NP_001137256.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07332
Protein Weight	93kDa
Shipping	Ice bag