

## Fgb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P02675

Item Number: CM0015607



### Introduction

Rabbit polyclonal antibody against human fibrinogen beta chain (FGB) validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for hemostasis, thrombosis, and wound healing research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | FGB Rabbit pAb                                                                                                    |
| Remarks/Alias    | HEL-S-78p; FGB                                                                                                    |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 2244                                                                                                              |
| GeneID           | 2244                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 111-222 of human FGB (NP_005132.2). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P02675                                                                                                            |
| Protein Weight   | 56kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |