

Fgd4 Rabbit Pab (Polyclonal Antibody) For Wb, Elisa - Q96M96

Item Number: CM0030582

Introduction

FGD4 Rabbit pAb is a polyclonal antibody targeting human FGD4, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for cytoskeleton dynamics and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	FGD4 Rabbit pAb
Remarks/Alias	CMT4H; FRABP; ZFYVE6; FGD4
Species	Human
GeneID (Human)	121512
GeneID	121512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-198 of human FGD4 (NP_640334.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96M96
Protein Weight	87kDa
Shipping	Ice bag