

Fgf14 Rabbit Pab For Wb, Elisa - Q92915

Item Number: CM0030704

Introduction

Rabbit polyclonal antibody against human fibroblast growth factor 14 (FGF14), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying nervous system development and function.

[Learn More](#)

Parameter	Value
Product Name	FGF14 Rabbit pAb
Remarks/Alias	FHF4; NYS4; FHF-4; SCA27; FGF-14; SCA27A; SCA27B; FGF14
Species	Human
Gene ID (Human)	2259
Gene ID	2259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human FGF14 (NP_787125.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92915
Protein Weight	28kDa
Shipping	Ice bag