

Fgf2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09038

Item Number: CM0003442

Introduction

FGF2 Rabbit Monoclonal Antibody validated in WB, IHC-P, IF/ICC, ELISA. Detects human FGF2 (bFGF), a heparin-binding growth factor involved in cell proliferation, differentiation, and angiogenesis. Ideal for oncology and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	FGF2 Rabbit mAb
Remarks/Alias	BFGF; FGFB; FGF-2; HBGF-2; F2
Species	Human
GeneID	2247
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human FGF2 (NP_001997.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P09038
Protein Weight	31kDa
Shipping	Ice bag