

Fgf23 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Gzv9

Item Number: CM0027147



Introduction

FGF23 Rabbit pAb validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Regulates phosphate and vitamin D metabolism; essential for bone disorder research. 28kDa target.

[Learn More](#)

Parameter	Value
Product Name	FGF23 Rabbit pAb
Remarks/Alias	ADHR; FGFN; HYPF; HFTC2; HPDR2; PHPTC; FGF23
Species	Human
Gene ID (Human)	8074
Gene ID	8074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-251 of human FGF23 (NP_065689.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV9
Protein Weight	28kDa
Shipping	Ice bag