


Fgfr2 Rabbit Polyclonal Antibody For Ihc-P And Elisa - P21802

Item Number: CM0022587



Introduction

Rabbit polyclonal antibody targeting human FGFR2, suitable for IHC-P and ELISA. Cross-reacts with mouse and rat. Ideal for developmental biology, oncology, and skeletal disorder research.

[Learn More](#)

Parameter	Value
Product Name	FGFR2 Rabbit pAb
Remarks/Alias	BEK; JWS; BBDS; CEK3; CFD1; ECT1; KGFR; TK14; TK25; BFR-1; CD332; K-SAM; FGFR2
Species	Human
Gene ID (Human)	2263
Gene ID	2263
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21802
Protein Weight	92kDa
Shipping	Ice bag