

Fhl1 Rabbit Pab For Wb And Elisa - Q13642

Item Number: CM0024864



Introduction

This rabbit polyclonal antibody recognizes the four and a half LIM domains protein 1 (FHL1/SLIM), a muscle development regulator. Validated in WB and ELISA, it cross-reacts with mouse, making it suitable for skeletal muscle and hypertrophy studies.

[Learn More](#)

Parameter	Value
Product Name	FHL1 Rabbit pAb
Remarks/Alias	KYOT; SLIM; FCMSU; FHL-1; FHL1A; FHL1B; FLH1A; SLIM1; XMPMA; RBMX1A; RBMX1B; SLIM-1; SLIMMER; FHL1
Species	Human
GeneID	2273
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-296 of human FHL1 (NP_001153171.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13642
Protein Weight	36kDa
Shipping	Ice bag