


Fip200 Rabbit Pab - Q8Tdy2

Item Number: CM0015282



Introduction

FIP200 Rabbit pAb is a polyclonal antibody detecting human, mouse, and rat RB1CC1. Validated in WB, IHC-P, ELISA; ideal for autophagy and cancer research.

[Learn More](#)

Parameter	Value
Product Name	FIP200 Rabbit pAb
Remarks/Alias	CC1; ATG17; FIP200; PPP1R131
Species	Human
GeneID (Human)	9821
GeneID	9821
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1280-1400 of human FIP200 (NP_055596.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDY2
Protein Weight	183kDa
Shipping	Ice bag