

Fitc Anti-Chicken Cd4 Monoclonal Antibody For Flow Cytometry - Q9W6V7

Item Number: CM0004644

Introduction

FITC-conjugated rabbit monoclonal antibody targeting chicken CD4 (55 kDa). Suitable for flow cytometry on chicken lymphocytes. Mapped to the extracellular domain of the T-cell surface protein.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Chicken CD4 mAb
Remarks/Alias	CD4
Species	Chicken
Gene ID (Human)	920
Gene ID	395362
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-423 of chicken CD4 (NP_989980.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Chicken
SWISS	Q9W6V7
Protein Weight	55 kDa
Shipping	Ice bag