

## Fitc Anti-Mouse Cd3E Monoclonal Antibody - Cd3E

Item Number: CM0003976

### Introduction

FITC-conjugated monoclonal antibody specifically recognizing mouse CD3e (CD3 epsilon) for flow cytometry applications. Sourced from Armenian hamster, this reagent enables precise detection of T lymphocytes in mouse samples.

[Learn More](#)

| Parameter        | Value                             |
|------------------|-----------------------------------|
| Product Name     | FITC Anti-Mouse CD3e mAb          |
| Remarks/Alias    | CD3; T3e; CD3epsilon              |
| Species          | Mouse                             |
| Gene ID          | 12501                             |
| Immunogen        | Recombinant protein of Mouse CD3e |
| Source           | Armenian Hamster                  |
| Category         | Monoclonal Antibodies             |
| Application      | FC                                |
| Cross Reactivity | Mouse                             |
| Shipping         | Ice bag                           |