

Fitc Rabbit Anti-Goat Igg (Fc) Fragment Monoclonal Antibody - Goat Igg Fc

Item Number: CM0009886

Introduction

High-specificity FITC-conjugated rabbit monoclonal antibody targeting goat IgG Fc fragment, validated for flow cytometry. Ideal for detecting goat IgG in research and diagnostic immunoassays.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit Anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein Recombinant protein of goat IgG (Fc) fragment.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag