


Fitc Rabbit Anti-Human Cd2 Monoclonal Antibody For Fc - P06729

Item Number: CM0008149



Introduction

FITC-conjugated rabbit anti-human CD2 monoclonal antibody for flow cytometry (FC). Specific for human CD2 (T11/SRBC/LFA-2), a 39kDa T-cell surface antigen. Validated for FC. Ideal for T-cell and NK cell immunophenotyping and adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Human CD2 mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
GeneID (Human)	914
GeneID	914
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag