


Fitc Rabbit Anti-Human Cd44 Monoclonal Antibody For Fc - P16070

Item Number: CM0008861



Introduction

FITC-conjugated rabbit monoclonal antibody targeting human CD44 for flow cytometry. Validated for human samples. Ideal for studying cell adhesion, migration, and immune homing. SWISS: P16070.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Human CD44 mAb
Remarks/Alias	CD44; CDW44; CSPG8; ECMR-III; HCELL; HUTCH-I; IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CD44 antigen
Species	Human
GeneID (human)	960
GeneID	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-178 of human CD44 (NP_000601.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag