

Fitc Rabbit Anti-Human Cd7 Monoclonal Antibody For Flow Cytometry - P09564

Item Number: CM0006111

Introduction

FITC-conjugated rabbit monoclonal antibody targeting human CD7 (GP40, Tp40, LEU-9). Designed for flow cytometric detection of T-cell and NK cell populations in immunophenotyping and diagnostic workflows.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Human CD7 mAb
Remarks/Alias	GP40; TP41; Tp40; LEU-9
Species	Human
GeneID (Human)	924
GeneID	924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag