


Fitc Rabbit Anti-Human Cd99 Monoclonal Antibody For Fc - P14209

Item Number: CM0007864



Introduction

FITC-conjugated rabbit monoclonal antibody against human CD99 (MIC2) for flow cytometry. Targets the 19 kDa CD99 antigen involved in T-cell adhesion, leukocyte extravasation, and rosette formation. Cross-reacts with human CD99.

[Learn More](#)

Parameter	Value
Product Name	FITC Rabbit anti-Human CD99 mAb
Remarks/Alias	MIC2; HBA71; MIC2X; MIC2Y; MSK5X
Species	Human
GeneID	4267
Immunogen	Synthetic peptide corresponding to amino acids 24-123 of human CD99 (NP_002405.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14209
Protein Weight	19kDa
Shipping	Ice bag